Istruzioni d'uso

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IFLEX Form**THERMODUR**

Technical information

Clear and hard disc, easy to shape and finish, made with Polyethylenterephthalat-Glycol Copolyester (PET-G). For the use of this material, we recommend that you use one of the models of 1 FLEX FORM thermoformers.

PRODUCT FOR PROFESSIONAL USE

Applications

- Hard byte for bleaching.
- Hard byte for temporary matrices.
- Hard byte for temporary crowns and bridges.
- Hard byte for surgical spoons.
- Hard byte for retainers.

Thermoforming

- Build a good plaster model using hard plaster or good quality 3D resin.
- The plaster model should not be very high and therefore should be thinned as much as possible. The base of the model must have a smooth surface. Delete the parts of the model (vestibule, buccal floor) that can prevent proper thermoforming. Carefully remove the undercut areas and round off the sharp edges of the plaster.
- The resin model, printed with a 3D printer, it is recommended to make it as low as possible and in the shape of a U-shaped "horseshoe".
- During thermoforming there is always a thinning of the material compared to the original thickness. In general, for every cm of height of the model there is a loss of thickness of about 20%. For this reason it is necessary to lower the height of the model as much as possible or insert the model into the well of the thermoformer and cover it with stainless steel granulate, leaving only the area to be formatted uncovered.
- Heat the disc respecting the temperatures and times listed in the section "Temperatures and times"
- Lower the clamp plate of the 1 FLEX FORM thermoformer until the heated disc completely surrounds the plaster or 3D resin model.

- Activate the vacuum and in a few seconds you will get your byte.
- Cool the byte just realized before removing it from the model.
- To facilitate removal, cut the plastic material as close as possible to the edge of the model using scissors and cutters.
- While the patient is using the byte, it is recommended to sanitize it daily.

NOTE: The disc is equipped with a double thermoformable protective film which allows, if necessary, to maintain the original brilliance and luster of the byte after thermoforming.

Thickness and format

This disc is available in the following thicknesses and formats:

SQUARE

125 X 125 mm

- 0000610 - Thickness 0,5 mm - 10 pcs
- 0000611 - Thickness 0,8 mm - 10 pcs
- 0000612 - Thickness 1 mm - 10 pcs
- 0000613 - Thickness 1,5 mm - 10 pcs
- 0000614 - Thickness 2 mm - 5 pcs
- 0000615 - Thickness 3 mm - 5 pcs

ROUND

Ø 120 mm

- 0000610R - Thickness 0,5 mm - 10 pcs
- 0000611R - Thickness 0,8 mm - 10 pcs
- 0000612R - Thickness 1 mm - 10 pcs
- 0000613R - Thickness 1,5 mm - 10 pcs
- 0000614R - Thickness 2 mm - 5 pcs
- 0000615R - Thickness 3 mm - 5 pcs

ROUND

Ø 125 mm

- 0000610T - Thickness 0,5 mm - 10 pcs
- 0000611T - Thickness 0,8 mm - 10 pcs
- 0000612T - Thickness 1 mm - 10 pcs
- 0000613T - Thickness 1,5 mm - 10 pcs
- 0000614T - Thickness 2 mm - 5 pcs
- 0000615T - Thickness 3 mm - 5 pcs

Temperatures and times

Listed below are the temperature and the heating time to obtain an optimal thermoforming.

- Thickness 0,5 mm - **157 C°** or **80 seconds**
- Thickness 0,8 mm - **162 C°** or **89 seconds**
- Thickness 1 mm - **165 C°** or **110 seconds**
- Thickness 1,5 mm - **152 C°** or **140 seconds**
- Thickness 2 mm - **155 C°** or **142 seconds**
- Thickness 3 mm - **147 C°** or **145 seconds**

NOTE: Temperatures and heating times listed are approximate and have been calculated using the 1 FLEX FORM thermoforming. It is advisable to adjust the temperatures and times to your needs.

Symbology

Consult instruction for use

Medical device

Keep away from sunlight

Protect from rain and humidity

Warnings

In the event of a serious accident, contact: 1 FLEX TECHNOLOGY S.R.L. - Phone: +39 0775 202311 - E-mail: info@1flex.it

BART MEDICAL S.R.L.

Distributed by

1 FLEX TECHNOLOGY S.R.L.

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03100 Frosinone (FR)

+39 0775 202311

info@1flex.it

www.1flex.it

Support

For further information and technical support please contact us by e-mail at info@1flex.it